Table of Contents

README	1.1
INSTALL	1.2
Users Guide	1.3
firmware	1.3.1
Decode	1.3.2
GUI	1.3.3
OnlineStattics	1.3.4
MakeEvent	1.3.5
FrontPanel	1.3.6
Programmers Guide	1.4
XIA API	1.4.1

README

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Base Setup

The Base Setup dialog box contains a table with 16 channels (ch # 0 to 15). Each channel has the following parameters: DCOffset, BLcut, Baseline[%], TraceDelay[us], TraceLength[us], Gain, Sign, GC, TC, EQS, ECT, and ERB. The values for DCOffset, BLcut, Baseline[%], TraceDelay[us], and TraceLength[us] are consistent across all channels. The Gain is set to 'Larger' and Sign to '-' for all channels. The checkboxes for GC, TC, EQS, ECT, and ERB are all checked. Below the table, there is a Status field showing 'Ready', a Module # dropdown set to 0, and buttons for Load, Apply, and Cancel. At the bottom, there is a Select channel # dropdown set to 0, and buttons for Copy, AdjustOffset, and BLcutFinder.

ch #	DCOffset	BLcut	Baseline[%]	TraceDelay[us]	TraceLength[us]	Gain	Sign	GC	TC	EQS	ECT	ERB
0	0.611	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
1	0.608	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
2	0.611	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
3	0.608	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
4	0.620	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
5	0.611	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
6	0.625	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
7	0.617	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
8	0.611	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
9	0.614	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
10	0.608	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
11	0.593	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
12	0.599	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
13	0.593	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
14	0.608	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
15	0.582	0	10	1.000	3.000	Larger	-	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Status: Ready
Module # 0 Load Apply Cancel
Select channel # 0 Copy AdjustOffset BLcutFinder

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- 16 channel gate input
- 1 module gate input
- 1 not use current

J151-J155

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single-ended TTL external input signals

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XIA API

```
PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16InitSystem (
    unsigned short NumModules,    // total number of Pixie16 modules in the system
    unsigned short *PXISlotMap,   // an array containing the PXI slot number for ea
ch pixie16 module
    unsigned short OfflineMode ); // specify if the system is in offline mode

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ExitSystem (
    unsigned short ModNum );      // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadModuleInfo (
    unsigned short ModNum,        // module number
    unsigned short *ModRev,       // returned module revision
    unsigned int *ModSerNum,      // returned module serial number
    unsigned short *ModADCBits,   // returned module ADC bits
    unsigned short *ModADCMSPS ); // returned module ADC sampling rate

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16BootModule (
    char *ComFPGAConfigFile,      // name of communications FPGA configuration file
    char *SPFPGAConfigFile,      // name of signal processing FPGA configuration f
ile
    char *TrigFPGAConfigFile,     // name of trigger FPGA configuration file
    char *DSPCodeFile,           // name of executable code file for digital signa
l processor (DSP)
    char *DSPParFile,            // name of DSP parameter file
    char *DSPVarFile,            // name of DSP variable names file
    unsigned short ModNum,        // pixie module number
    unsigned short BootPattern ); // boot pattern bit mask

// Acquire ADC traces in single or multiple modules
// Use this function to acquire ADC traces from Pixie-16 modules. Specify the modul
e using ModNum. If ModNum is set to be less than the total number of modules in the
system, only the module specified by ModNum will have its ADC traces acquired. But
if ModNum is equal to the total number of modules in the system, then all modules
in the system will have their ADC traces acquired.
// After the successful return of this function, the DSP's internal memory will be
filled with ADC trace data. A user's application software should then call another
function Pixie16ReadSglChanADCTrace to read the ADC trace data out to the host comp
uter, channel by channel.
PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16AcquireADCTrace (
    unsigned short ModNum );      // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadSglChanADCTrace (
    unsigned short *Trace_Buffer, // trace data
    unsigned int Trace_Length,    // trace length
    unsigned short ModNum,        // module number
    unsigned short ChanNum );     // channel number
```

```
PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16IMbufferIO (
    unsigned int *Buffer,          // buffer data
    unsigned int NumWords,         // number of buffer words to read or write
    unsigned int Address,          // buffer address
    unsigned short Direction,      // I/O direction
    unsigned short ModNum );       // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16EMbufferIO (
    unsigned int *Buffer,          // buffer data
    unsigned int NumWords,         // number of buffer words to read or write
    unsigned int Address,          // buffer address
    unsigned short Direction,      // I/O direction
    unsigned short ModNum );       // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16StartListModeRun (
    unsigned short ModNum,         // module number
    unsigned short RunType,        // run type
    unsigned short mode );         // run mode

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16StartHistogramRun (
    unsigned short ModNum,         // module number
    unsigned short mode );         // run mode

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16CheckRunStatus (
    unsigned short ModNum );       // Pixie module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16EndRun (
    unsigned short ModNum );       // Pixie module number

PIXIE16APP_EXPORT double PIXIE16APP_API Pixie16ComputeInputCountRate (
    unsigned int *Statistics,
    unsigned short ModNum,
    unsigned short ChanNum );

PIXIE16APP_EXPORT double PIXIE16APP_API Pixie16ComputeOutputCountRate (
    unsigned int *Statistics,
    unsigned short ModNum,
    unsigned short ChanNum );

PIXIE16APP_EXPORT double PIXIE16APP_API Pixie16ComputeLiveTime (
    unsigned int *Statistics,
    unsigned short ModNum,
    unsigned short ChanNum );

PIXIE16APP_EXPORT double PIXIE16APP_API Pixie16ComputeProcessedEvents (
    unsigned int *Statistics,
    unsigned short ModNum );

PIXIE16APP_EXPORT double PIXIE16APP_API Pixie16ComputeRealTime (
    unsigned int *Statistics,
```

```
    unsigned short ModNum );

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16complexFFT (
    double *data,
    unsigned int length );

PIXIE16APP_EXPORT unsigned short PIXIE16APP_API APP16_TstBit (
    unsigned short bit,
    unsigned short value );

PIXIE16APP_EXPORT unsigned short PIXIE16APP_API APP16_SetBit (
    unsigned short bit,
    unsigned short value );

PIXIE16APP_EXPORT unsigned short PIXIE16APP_API APP16_ClrBit (
    unsigned short bit,
    unsigned short value );

PIXIE16APP_EXPORT unsigned int PIXIE16APP_API APP32_SetBit (
    unsigned short bit,
    unsigned int value );

PIXIE16APP_EXPORT unsigned int PIXIE16APP_API APP32_ClrBit (
    unsigned short bit,
    unsigned int value );

PIXIE16APP_EXPORT unsigned int PIXIE16APP_API APP32_TstBit (
    unsigned short bit,
    unsigned int value );

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16SetDACs (
    unsigned short ModNum );

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ProgramFippi (
    unsigned short ModNum );

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16AdjustOffsets (
    unsigned short ModNum );

// Acquire baselines from a module
PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16AcquireBaselines (
    unsigned short ModNum );    // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadSglChanBaselines (
    double *Baselines,          // returned baselines values
    double *TimeStamps,         // time stamp for each baseline value
    unsigned short NumBases,     // number of baseline values to read
    unsigned short ModNum,       // module number
    unsigned short ChanNum );    // channel number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16RampOffsetDACs (
```

```

    double *DCValues,          // returned DC offset values
    unsigned short NumDCVals,   // number of DC values to read
    unsigned short ModNum );

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ControlTaskRun (
    unsigned short ModNum,      // Pixie module number
    unsigned short ControlTask, // Control task number
    unsigned int   Max_Poll );  // Timeout control in unit of ms for control task
run

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16BLcutFinder (
    unsigned short ModNum,      // Pixie module number
    unsigned short ChanNum,     // Pixie channel number
    unsigned int   *BLcut );    // BLcut return value

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16TauFinder (
    unsigned short ModNum,      // Pixie module number
    double         *Tau );      // 16 returned Tau values, in 祢

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16WriteSglModPar (
    char *ModParName,          // the name of the module parameter
    unsigned int   ModParData, // the module parameter value to be written to th
e module
    unsigned short ModNum );    // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadSglModPar (
    char *ModParName,          // the name of the module parameter
    unsigned int   *ModParData, // the module parameter value to be read from the
module
    unsigned short ModNum );    // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16WriteSglChanPar (
    char *ChanParName,         // the name of the channel parameter
    double ChanParData,        // the channel parameter value to be written to t
he module
    unsigned short ModNum,      // module number
    unsigned short ChanNum );   // channel number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadSglChanPar (
    char *ChanParName,         // the name of the channel parameter
    double *ChanParData,       // the channel parameter value to be read from th
e module
    unsigned short ModNum,      // module number
    unsigned short ChanNum );   // channel number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadHistogramFromModule (
    unsigned int   *Histogram,  // histogram data
    unsigned int   NumWords,    // number of words to be read out
    unsigned short ModNum,      // module number
    unsigned short ChanNum);    // channel number

```

```
PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadStatisticsFromModule (
    unsigned int    *Statistics,    // run statistics data
    unsigned short  ModNum );      // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16SaveHistogramToFile (
    char *FileName,                // histogram data file name
    unsigned short  ModNum);       // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16GetModuleEvents (
    char *FileName,                // the list mode data file name (with complete pa
th)
    unsigned int *ModuleEvents ); // receives number of events for each module

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16GetEventsInfo (
    char *FileName,                // the list mode data file name (with complet
e path)
    unsigned int    *EventInformation, // to hold event information
    unsigned short  ModuleNumber);     // the module whose events are to be retrieved

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadListModeTrace (
    char *FileName,                // list mode data file name
    unsigned short *Trace_Data,    // list mode trace data (16-bit words)
    unsigned short NumWords,       // number of 16-bit words to be read out
    unsigned int   FileLocation);  // the location of the trace in the file

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadHistogramFromFile (
    char *FileName,                // the histogram data file name (with complete pa
th)
    unsigned int *Histogram,       // histogram data
    unsigned int  NumWords,        // number of words to be read out
    unsigned short ModNum,         // module number
    unsigned short ChanNum);       // channel number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16SaveDSPPParametersToFile (
    char *FileName );              // the DSP parameters file name (with complete pa
th)

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16LoadDSPPParametersFromFile (
    char *FileName );              // the DSP parameters file name (with complete pa
th)

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16CopyDSPPParameters (
    unsigned short BitMask,        // copy items bit mask
    unsigned short SourceModule,   // source module
    unsigned short SourceChannel,  // source channel
    unsigned short *DestinationMask ); // the destination module and channel bit m
ask

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadMSGFile (
    char *ReturnMsgStr );
```

```

PIXIE16APP_EXPORT unsigned int PIXIE16APP_API Decimal2IEEEFloating(double DecimalNumber);

PIXIE16APP_EXPORT double PIXIE16APP_API IEEEFloating2Decimal(unsigned int IEEEFloat
ingNumber);

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16SaveExternalFIFODataToFile (
    char *FileName,           // list mode data file name
    unsigned int *nFIFOWords, // number of words read from external FIFO
    unsigned short ModNum,    // module number
    unsigned short EndOfRunRead); // indicator whether this is the end of run read

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16RegisterIO (
    unsigned short ModNum,    // the Pixie module to communicate to
    unsigned int address,     // register address
    unsigned short direction, // either MOD_READ or MOD_WRITE
    unsigned int *value );    // holds or receives the data

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadCSR (
    unsigned short ModNum,
    unsigned int *CSR );

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16WriteCSR (
    unsigned short ModNum,
    unsigned int CSR );

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16CheckExternalFIFOStatus (
    unsigned int *nFIFOWords,
    unsigned short ModNum );

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ReadDataFromExternalFIFO (
    unsigned int *ExtFIFO_Data, // To receive the external FIFO data
    unsigned int nFIFOWords,    // number of words to read from external FIFO
    unsigned short ModNum );    // module number

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ComputeFastFiltersOffline (
    char *FileName,           // the list mode data file name (with complete path)
    unsigned short ModuleNumber, // the module whose events are to be analyzed
    unsigned short ChannelNumber, // the channel whose events are to be analyzed
    unsigned int FileLocation,    // the location of the trace in the file
    unsigned short RcdTraceLength, // recorded trace length
    unsigned short *RcdTrace,     // recorded trace
    double *fastfilter,           // fast filter response
    double *cfd );               // cfd response

PIXIE16APP_EXPORT int PIXIE16APP_API Pixie16ComputeSlowFiltersOffline (
    char *FileName,           // the list mode data file name (with complete path)

```

```

te path)
    unsigned short ModuleNumber,          // the module whose events are to be analyzed

    unsigned short ChannelNumber,         // the channel whose events are to be analyzed
ed
    unsigned int   FileLocation,           // the location of the trace in the file
    unsigned short RcdTraceLength,        // recorded trace length
    unsigned short *RcdTrace,             // recorded trace
    double         *slowfilter );         // slow filter response

// Add by Hongyi Wu
PIXIE16APP_EXPORT int PIXIE16APP_API HongyiWuPixie16ComputeSlowFiltersOffline (
    char           *FileName,              // the list mode data file name (with complete path)
    unsigned short ModuleNumber,           // the module whose events are to be analyzed

    unsigned short ChannelNumber,         // the channel whose events are to be analyzed
ed
    unsigned int   FileLocation,           // the location of the trace in the file
    unsigned short RcdTraceLength,        // recorded trace length
    unsigned short *RcdTrace,             // recorded trace
    double         *slowfilter,           // slow filter response
    unsigned int   bl,
    double         sl,
    double         sg,
    double         tau,
    int            sfr,
    int            pointtobl );

```